

Standard of Care – Bloods to Idexx

1. Blood sample collection kit is located in both consult rooms, and prep room. Consists of
 - x Red 2ml clotting
 - x Purple 2ml Edta – invert once blood is in tube
 - 1 x glass slide
2. Once blood is collected and tubes are filled up to correct amount (see black line on tubes) print ID label from Vetlink and label blood tubes
3. Fill in Idexx blood form located in prep room; top RH cupboard
4. Call Idexx **1300 443399** press **1**, then prompt will notify on confirmation
5. Steve/reception will print out con note from email
 - One con note to go on parcel sleeve (toll bag)
 - Put animal name and last name on con note and place in 2nd drawer under printer out the front
6. Have all blood tubes labelled and in bio-hazard bag
 - Get white container located in bloods room, bottom R corner cupboard
 - Get 2 small ice packs, wrap in newspaper, then fold bloods up so all can fit into container
 - Get Toll bag, clear sleeve and con note
7. Put tick on Fridge (to alert that there are bloods to go)
8. **MAKE SURE EVERYTHING IS READY BY 2.30 FOR TOLL COLLECTION**